

FLUID CONTROL SOLIDIFIERS

For potentially infectious
blood and bodily fluid waste



Red Z® and Green-Z® Fluid Control Solidifiers

Healthcare facilities must focus on maintaining high levels of worker safety while reducing operating costs. Help improve surgical room turnover and worker safety when handling or transporting liquid medical waste with the use of Red Z® and Green-Z® bodily fluid solidifiers.



Fast-Acting

Contains superabsorbent polymers (SAPs)



Shipping-Friendly

Helps comply with DOT regulations



No Splashing

Eliminates splashing and aerosolization



Doesn't Expire*

Shelf-stable in a cool, dry place

Industries at Risk in Handling Fluid Waste

- Healthcare Facilities
- Research Centers
- Dentists & Orthodontists
- Mortuaries
- Autopsy Centers
- Testing Laboratories
- Blood Banks
- Collection Services
- Nursing Homes
- Tattoo Parlors
- Correctional Facilities
- Police & Fire Departments

**See following pages for specific
product information**


*Products are shelf-stable and do not require an expiration date. For the product to remain shelf-stable, it must be kept in a cool, dry place and in its original packaging. If contents are hard/compact, shake the product to loosen before use.



Safetec of America, Inc.
1-800-456-7077 | safetec.com

Solidifiers are Available in Three Packaging Options:

- ① Single & Multi-Use Bottles
- ② Single-Use Pouches
- ③ Drop-In Pacs



“Simple waste management measures, such as effective confinement of waste and safe handling, can already dramatically reduce health risks...” – World Health Organization

① Single & Multi-Use Bottles

- Pour into porthole of suction canisters
- Also great for urology bags and kick buckets
- Available in a variety of sizes to solidify up to 22,000cc of liquid medical waste

Item#	Description	Case Qty.
41000	Red Z [®] Up to 1,500cc Bottle	75/case
41005	Red Z [®] Up to 3,000cc Bottle	75/case
41105	Red Z [®] Up to 5,000cc Bottle	24/case
41107	Red Z [®] Up to 22,000cc Bottle	12/case
42005	Green-Z [®] Up to 3,000cc Bottle	75/case
42003	Green-Z [®] Up to 22,000cc Bottle	12/case



Great for suction canisters, urology bags, and kick buckets

Directions for Use

Suction Canisters:

1. Release proper amount of solidifier through port opening of suction canister (before or after surgical procedure).
2. Do not allow product packaging to contact contaminated material.
3. Allow sufficient time for encapsulation before handling.

Biohazard Bags or Kick Buckets:

1. Add solidifier to the bottom of a biohazard bag or other container to control fluid before or after spills occur.
2. Clean and disinfect contaminated areas after use.
3. Dispose of infectious waste in accordance with all local, state and federal regulations.

② Pour-In Pouches

- Single-use, pre-measured powder
- Easy-pour angled pouch is perfect for narrow openings of suction canisters
- Snip with scissors and pour in
- Red Z® contains chlorine to deodorize and control odor



Item#	Description	Case Qty.
41131	Red Z® 1oz. Pour-In Pouch	100/case
41132	Red Z® 2oz. Pour-In Pouch	100/case

Directions for Use

Suction Canisters:

1. Use scissors to snip corner of pouch and release proper amount of solidifier through port opening of suction canister (before or after surgical procedure).
2. Discard product packaging.
3. Allow sufficient time for solidification and encapsulation of contents before handling.

Biohazard Bags or Kick Buckets:

1. Add solidifier to the bottom of a biohazard bag or other container to control fluid before or after spills occur.
2. Clean and disinfect contaminated areas after use.
3. Dispose of infectious waste in accordance with all local, state and federal regulations.

③ Drop-In Pacs

- Great for pre-loading and post-loading canisters
- Easy and ready-to-use pre-measured pouches
- Fast-acting water-soluble paper technology
- Contents are non-toxic

Item#	Description	Case Qty.
41108	Red Z® Drop-In (Up to 750cc)	300 per cs.
41175	Red Z® Drop-In (Up to 1500cc)	100 per cs.
42008	Green-Z® Drop-In (Up to 750cc)	300 per cs.
42175	Green-Z® Drop-In (Up to 1500cc)	100 per cs.



**Specifically designed to fit in
the porthole of a canister**

Directions for Use

Suction Canisters:

1. Insert entire drop-in pac through port opening of suction canister (before or after surgical procedure).
2. Drop-in Pac will rupture within a few minutes upon contact of bio-fluid.
3. Allow sufficient time for solidification and encapsulation of contents before handling.

Biohazard Bags or Kick Buckets:

1. Place drop-in pac in the bottom of a biohazard bag or other container to control fluid before or after spills occur.
2. Clean and disinfect contaminated areas after use. Dispose of infectious waste in accordance with all local, state and federal regulations.