# ALLOGEL® MAGENTA™





← SmartPhone Scan





Nitrile Examination Glove with AlloGel® Inner Coating, ULTRA SOFT-STRETCH™, Powder Free, Micro-Textured, Beaded Cuff, Non-Sterile, Magenta™ Color

#### MATERIAL SAFETY DATA SHEET (MSDS): HEALTH HAZARD DATA

The chemical formulation of the gloves and the chlorination process do not contain any substances known to be harmful to the user. You are required to follow any cautionary and warning labels on the product packaging. Some users may be allergic to components used in the manufacturing process. If this is the case, you must immediately stop using this product and consult your medical doctor. These gloves have undergone BIOCOMPATABILITY tests as recommended by the US FDA for Primary Skin Irritation and Dermal Sensitization and were found to be non-irritant and non-sensitive. For more information, please call us in the USA at 651-322-9112.

# **Technical Specifications**

**Storage**: Avoid direct sunlight. Store in a cool and dry place. EN-455

Tensile Strength (Mpa): Before aging Min. 14 Min. 6N

After aging at 70°C, 7 days Min. 14 Min. 6N

Ultimate Elongation %: Before aging Min. 500

After aging at 70°C, 7 days Min. 400

Size	XS	S	М	L	XL	XXL
Width (±5 mm)	70	80	95	110	115	125
Length (±5 mm)	240	240	240	240	240	240
Thickness						
Cuff (±0.03 mm)	0.07	0.07	0.07	0.07	0.07	0.07
Palm (±0.03 mm)	0.08	0.08	0.08	0.08	0.08	0.08
Finger (±0.05 mm)	0.12	0.12	0.12	0.12	0.12	0.12
Mil	3	3	3	3	3	3

# Standards & Certifications & Chemical Composition

ASTM D 6124 Residual powder in medical gloves EN 374
ASTM D 6319 Nitrile Exam Gloves for medical application EN 420
ASTM D 5151 Detection of holes in medical gloves EN 455
ASTM D 412 Tension test ISO 2859
ASTM F 1671 Viral Penetration

EN 455 EU Medical Gloves ISO 2859 Sampling Procedures

Complies with the Class 1 requir

Class 1
Complies with the Class 1 requirements of the European Medical Device Directive



EN374-2003



Complies with EN-374-1 & 2, Level 2 European Medical Device Directive The Chemical Composition of this glove includes some or all of the following ingredients:

Protection from Chemicals & Micro-organisms

**General Requirements for Gloves** 

Acrylonitrile Butadiene, Potassium Hydroxide, Sulphur, ZDBC, ZDEC, Titanium Dioxide, Zinc Oxide, , Calcium Carbonate, Calcium Nitrate, Sodium Hydroxide, Lauric Acid, Nitric Acid, Sodium Hypochlorite, Ethoxylate fatty alcohol, NSFP sodium salt, Collagen, Allantoin, Silicone Emulsion. Magenta pigment.

XL

**ORDERING** 

Item#

ALLOGENA-XS

S ALLOGENA-S

M

**ALLOGENA-M** 

ALLOGENA-L

ALLOGENA-XL

XXL ALLOGENA-XXL

Packaging: XS-L: 100 gloves per box by weight; 10 boxes per case; XL-XXL: 90 gloves per box by weight; 10 boxes per case

## Distributed by:

AMERICAS
Preventive Care, Inc.
www.PerfectGloves.com
+1 651 322 9199

EUROPE EuroDispo
EuroDispo A/S
www.EuroDispo.eu
+45 7582 0303

eXAsia Sdn Bhd www.eXAsia.asia +60 5 808 1115

### PATENTS & TRADEMARKS:

This product may be covered by one or more pending or issued patients, including U.S. Patent 8365314, European Patent (FPO) FP 2 081 454 B1.

AlloGena<sup>m,</sup> AlloGel<sup>a</sup>, AlloGel<sup>a</sup>, I, AlloGel<sup>a</sup> FIT, AlloGel<sup>a</sup> FIT, I, CLEAN-PERFECT<sup>m</sup>, CollaGel<sup>a</sup>, EXAM-PERFECT<sup>a</sup>, KLINI-PERFECT<sup>m</sup>, MICRO-PERFECT<sup>m</sup>, NITRI-PERFECT<sup>m</sup>, NITRID-PERFECT<sup>a</sup>, SAWCARE<sup>a</sup>, SYNTHO-PERFECT<sup>a</sup>, The PERFECT<sup>a</sup>, AND ALLO PERFECT<sup>a</sup>, SAWCARE<sup>a</sup>, SYNTHO-PERFECT<sup>a</sup>, AND ALLO PERFECT<sup>a</sup>, SAWCARE<sup>a</sup>, SYNTHO-PERFECT<sup>a</sup>, AND ALLO PERFECT<sup>a</sup>, SAWCARE<sup>a</sup>, Visible Wave<sup>m</sup>, www. AlloGena.com<sup>m</sup>, the Barcode design, the box design, the color COPPER<sup>a</sup> and the color MAGENIA<sup>m</sup> are trademarks of their respective owners. \*\* please note: patents and registered trademarks are specific to regions, countries, states or other defined areas. A specific patent or registered mark is always only assigned to cover a specific limited geographic area. If you wish to obtain the coverage of a mark, please contact ip@xelacorp.com & request the coverage details for the specific mark(s). © 2013 Preventive Care, Inc.